AUTONOM HEALTH

# SELF MONITORING for HEALTH & SPORTS

HRV QUICKGUIDE

# MY AUTONOM HEALTH & SPORTS

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

# MY AUTONOM HEALTH & SPORTS

*Visualize your health!*

**Heart Rate Variability measurements are very simple and amazingly meaningful:**

With minimal effort, they give precise and objective information about your physical and mental state, health, lifestyle and performance.

ENGLISH

**Instructions for a HRV measurement:**  
assistance to perform HRV measurements  
with Bluetooth chest belts.

*Further information:*  
[www.autonomhealth.com](http://www.autonomhealth.com)

## CONTENT

- Step by step to your measurement 36-45
- Technical details 46-48
- Maintenance and care of the measuring system 49-51
- Background knowledge about HRV 52-55
- LifeFire® Gallery 56-57
- Support with questions and problems 58-59

## EVERYTHING AT A GLANCE

1

**REGISTRATION IN THE ANALYSIS PORTAL**

*<https://portal.autonomhealth.com>*

2

**REDEEM VOUCHER CODE AND ACTIVATE  
TARIFF**

3

**DOWNLOAD THE APP**

4

**PUT ON THE CHEST BELT AND START THE  
APP ON YOUR PHONE**

5

**WATCH YOUR HEALTH AND STUDY THE  
RESULTS**

**Hint:** step 1-3 only necessary for the first measurement!

## A. PREPARATION

1

### REGISTRATION IN THE ANALYSIS PORTAL

*Your first step towards health*

Register in the analysis portal  
(<https://portal.autonomhealth.com>)  
with a personal account.

**Hint:** If you are already registered with a personal account, skip step 1 and log in. A personal account can only be used by ONE person!

## 2

## REDEEM A VOUCHER CODE

*A valid tariff allows everything*

The voucher code for your software package can be found in the box on the right of the sensor. Enter this code into your account on the analysis portal (<https://portal.autonomhealth.com>) on the page “My Account/My Order”. The tariff is activated and usable immediately after the redemption of the voucher.

The voucher activates the usage of *My Autonom Health & Sports* and includes:

- graphic analyses, results, explanations and recommendations of the ten most significant health and lifestyle factors
- detailed recommendation based on your measurements
- your 5-minute accurate performance analysis during the measurement
- the entire sports analysis program for your training

3

## DOWNLOAD THE APP

*Your health coach on your phone*

Download the *My Autonom Health App* for free from the App or Play Store. Log in or register with your access data.

**Hint:** If you have already registered in the app before redeeming the tariff, please log out and then log in again to activate the tariff in the app. (low battery, etc.). In order to always work with the latest app version, we recommend setting the automatic update on the mobile phone.

© AUTONOM HEALTH 

## B. MEASUREMENT

1

### ATTACH THE CHEST BELT

*The measuring device, which makes your heart speak*

- Click the sensor (black, round middle piece with the Autonom Health logo) with the push buttons on the belt. The logo on the sensor should be readable from the front.
- Attach the chest belt with the electrodes facing inwards around your chest and underneath your breast.
- Adjust the belt length so that the belt does not slip and the electrodes are attached closely.

## STARTING APP AND MEASUREMENT


### *The health App as a self-explanatory support*

Activate the **My Autonom Health App** on your smartphone.

Before first use, the sensor must be connected to the App. Select your sensor in the list of compatible devices and connect it to the App.

Once the pulse rate (bpm) is displayed, the sensor is connected to the App and the measurement can be started.

**Hint:** The App will alert you immediately if there is a problem (insufficient connection with the chest strap, low battery, etc.). The coupling of the Bluetooth connection must take place directly in the app and not on the smartphone!

## ON THE WAY TO YOUR HEALTH

*Watch your LifeFire® grow*

If there is an active internet and Bluetooth connection, your measured data is constantly sent to the analysis portal and gets analyzed there periodically. That's how you can follow the process of your HRV measurement permanently both in the App and in the analysis portal and you can always know, if your performance is correct or should be changed.

**Push notification:** Intelligent Pop-Up hints assist during the measurement. They give impulses, warn against overload and encourage you to improve your lifestyle. These push notifications can be turned on/off on request in the settings of your phone.

**Important Hint:** Thanks to the internal memory in the sensor, a permanent Bluetooth connection to the smartphone is not required, this means Data from a longer measurement unit without smartphone are buffered in the sensor and read out later. Transmission and analysis of this stored data can take several minutes!

4


## FINISH YOUR MEASUREMENT

*All results, recommendations and explanations are already waiting in the analysis portal*

When the measurement process is finished press the „End measurement button“. Then you can retrieve your results in the analysis portal.

*(<https://portal.autonomhealth.com>)*

**Note:** Make sure that all data are read out before ending the measurement! Otherwise, data will be lost and cannot be restored.

## C. MEASUREMENT RESULTS

1

### Open your analytical report

*Your health in all details*

Many results are available directly in the App. In the analysis portal you will find several detailed evaluations as well as detailed explanations and individual recommendations.

Discover analysis details on the five major **health determinants** like health status, biological age, performance potential, stress management, and burnout resistance, as well as the five decisive **lifestyle factors**: sleep, recovery, mental and physical performance and nutritional habits, as well as the complete sports analysis program for your exercise control. The performance bar shows your vitality state precisely up to 5 minutes and specifically for all your activities.

## D. CHOSEN TARIFF EXPIRED

### How to proceed

You are enthusiastic about the opportunity to take charge of your own health and still want to find out what is good for you. Then you have two easy possibilities to order another tariff:

1

Order directly in your private account in the analysis portal (<https://portal.lebensfeuer.com/payscales/index/payscales>) on the page **“My Account/Order now”**.

OR

2

Order directly in the App. With Android, you will be forwarded to the analysis portal with a link and then you can order in your personal account.

## E. THE BEST TIMING FOR A MEASUREMENT

For the first measurement, a daily routine with a night in your own bed is best. You should stay in your usual professional and private environment, sleep in your own bed and live quite normally throughout the day. Then there are hardly any excuses for not immediately starting measurements ...

Once you have collected your “basic data”, it’s time to experiment. Since not everyone reacts to situations, activities, food, etc. in the same way, thanks to My Autonom Health, it is finally up to you to find out what is good for you or what you might be well advised to have less of in your everyday life.

# Technical details

## Information about the sensor

- 9-axis motion sensor: 3x accelerometer, 3x gyroscope, 3x magnetometer
- Pulse/heart rate
- Integrated data storage
- Bluetooth Smart compatible
- Waterproof up to 30 m / 100 ft

## Dimensions

- Diameter: 36.6 mm (1.44 inches)
- Thickness: 10.6 mm (0.42 inches)
- Weight: 10 g including battery

## MCU

- Nordic Semiconductor nRF52832 consists of 32-Bit ARM® Cortex®-M4 with 64 KB RAM and 512 KB FLASH on the chip
- Integrated Bluetooth Low Energy 4.0-transmitter with 2.4 GHz

## Storage

- External storage: 3 Mbit EEPROM

**Hint:** The storage capacity is approx. 10-12 hours.

## Sensors

- Accelerometer:  $\pm 2 / \pm 4 / \pm 8 / \pm 16g$  full scale, sampling frequency: 12.5/26/52/104/208 Hz
- Gyroscope:  $\pm 125 / \pm 245 / \pm 500 / \pm 1000 / \pm 2000$  dps full scale, sampling frequency: 12.5/26/52/104/208 Hz
- Magnetometer:  $\pm 4 / \pm 8 / \pm 12 / \pm 16$  Gauss full scale value
- Temperature, precision:  $\pm 0,5^{\circ}C$  (max) from  $0^{\circ}C$  to  $+ 65^{\circ}C$  and  $\pm 1,0^{\circ}C$  (max) from  $- 40^{\circ}C$  to  $+ 125^{\circ}C$
- Heart frequency and single channel-ECG
- RR-interval sampling frequency 1.028 Hz

## Input/Output

- LED (red) at the front
- wire interface master (access over movesense-connector)
- Debug from UART, Debugging I / F, power und mass (in battery case)

## Energy source

- CR 2025 3 V lithium-coin cell battery

**Hint:** A quality battery lasts for approx. 6-8 long-term measurements. The app provides information when the battery power is low.

## Admission and conformities

- CE, FCC, IC, C-Tick, CMIIT
- Corresponds to REACH, RoHS
- Bluetooth 4.0

# Maintenance and care of the measuring system

## Cleaning

To clean the chest belt, use a soft cloth or, in the case of coarse impurities, a mixture of soap and water. Never use chemical cleaning agents, solvents, alcohol or abrasive materials.

**Attention:** The textile belt should be washed in the washing machine at 30-40°C from time to time. Any grease films or soiling can have a massive impact on the conductivity of the electrodes. Please do not forget to remove the sensor before washing as the sensor is not suitable for washing machines.

## Storage

After each use, please disconnect the sensor from the belt and store the two parts in a dry place. In this way, the battery should be more durable. Do not expose the sensor to extreme temperature changes over a long period.

### **Battery change**

Lift the battery compartment on the back of the sensor with a suitable object. Remove the empty battery (please dispose of it properly) and insert a new 3 V battery (type CR2025). Make sure that the positive pole (+) points upwards.

Make sure the cover is properly closed again (mind the notches) to prevent moisture from entering the sensor.

### **Guarantee**

The guarantee period for the measuring system, if properly used, is 12 months from the date of purchase (exclusively with proof of purchase) and is limited to material and manufacturing defects.

The guarantee expires in case of accidental damage as well as incorrect use of the product or improper battery replacement.

Excluded from the guarantee are batteries, normal wear and tear, damage caused by improper handling, accident, fall, impact, etc. Excluded from the warranty is the elastic electrode part.

**Attention:** The guarantee does not cover direct, indirect or consequential damages, costs, expenses or loss associated with the product.


### Heart Rate Variability - every heartbeat counts

Our heart is a fascinating muscle. It reacts directly to everything we experience externally, and even think and feel inside. The heartbeat is influenced by our inner clock, our breathing, our emotions and external influences. In healthy humans, the heart reacts continuously, like a high-tech instrument, to these external and internal signals with finely tuned changes (variations) in the heartbeat sequence. Hence the term “Heart Rate Variability“, abbreviated “HRV“.

#### **The astonishing:**

*My Autonom Health & Sports* measures your individual HRV. The result describes the ability of your heartbeat to adapt to different situations, the time interval from one beat to another.


This is a precise measure of the general adaptability of your organism and thus a measure of health, performance, regenerative ability and much more.

## Visualize your health


### What 10 g of technology and a picture of flames can do for your health and vitality

*My Autonom Health & Sports* acting as the information point of your heart clearly indicates your individual stress- and recovery-phases. Your heart adapts to the ever-changing external and internal conditions by shortening or lengthening the time between each heartbeat as needed. We use this physiological regularity to measure and understand health, fitness and lifestyle. To illustrate this dynamic insight into your health, a flame picture is created after the measurement.


## How does your flame picture come about?


The spectrogram, calculated from many thousands of data, resembles a fire: the more intense and denser the colors (blue, red, yellow, etc.) and the higher the “flames” blaze in the spectrogram, the more vital, younger and fit you are. If the diagram shows only small flames and little color contrasts, the picture indicates reduced vitality, little parasympathetic activity (regeneration, economy), exhaustion or even illness. The LifeFire® is a pictorial representation of your autonomic nervous system. Whether it’s a firework of vitality or an advanced burnout syndrome: everything is objectively documented because your body cannot lie.


## No two images are alike


*excellent sleep*


*bad sleep*

*burnout*

*top athlete*


*very good general state of health*

*fatigue*

*chronic alcohol misuse*

*good ability to concentrate*

## Support with questions & problems

### Questions? We can help!

#### Video clips

Immediate support can be found directly in the analysis portal! The **explanatory video** clips already answer many of your questions about your HRV measurements.

#### SUPPORT Button in the analysis portal

The **SUPPORT Button** in the green header of the analysis portal facilitates direct contact with our **support team** if problems or questions occur.

## Help from the Autonom Health® Community

Address your questions directly to our **expert forum** in the Autonom Health® Community. There you will receive free help and answers from our numerous **professionals** and **partners** of various health professions.

## Service Hotline


Got an urgent question?  
We are available by phone from  
**Monday to Friday from 9:00 to 17:00**  
under **+43 699 122 000 01.**

## Mail to the Service Center

You can also reach the service center via our email address **info@autonomhealth.com.**

Please expect a processing time of about 2 days for our answer. For specific questions about individual measurements, a fee will be charged. It covers the costs related to the time spent by our experts to provide the individual answer.

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